594 mm

49,5 mm

200 mm

TR

TEKNOMEK

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HYDROXYPROPYL METHYLLULOSE (VISCOTEK) INSTRUCTIONS FOR USE

Read the instructions for use cautiously!

1.PRODUCT INFORMATION

Product Name: Hydroxypropyl Methylcellulose

Product Model: See the catalog or web-site for the models

BUDI-DI : 8680548990226NM

Product Concentration Percentages: 2%, 2.4%


Size: 1 or 2 milliliters (in a glass syringe)

Shelf Life: Shelf life is two years. Stated on the product.

Class: IIb

Chemical Properties	Viscotek %2	Viscotek % 2.4
Methylcellulose Concentration Percentage	%2	%2.4
Molecular Weight	900.000-1.100.000 Da	900.000-1.100.000 Da
pH	6.8-7.6	6.8-7.6
Density	0.987-0.995 g/cm³	0.992-0.998 g/cm³
Concentration	20mg/ml	24mg/ml
Viscosity	10.000-20.000mPa.s	20.000-30.000mPa.s
Osmolarity	200-400mOsm/kg	200-400mOsm/kg
Endotoxicity	<0.2EU/ml	<0.2EU/ml
Refractive Index	1.337-1.338	1.337-1.338

Rheological Profile

Methylcellulose is a sterile, single use, surgically invasive viscoelastic solution. Methylcellulose is a clear, colorless, odorless and sterile solution for intraocular use. Does not contain any preservatives. It is indicated for use as an assistant material to ophthalmic surgery in any anterior segment procedures like cataract extraction and intraocular lens implantation. The solution contains hydroxypropyl methylcellulose in physiological saline solution at different concentrations.

2.PACKAGING and STERILIZATION

Methylcellulose is offered in 1 ml or 2 ml sterile single use glass syringes. The product is offered with a 23-27 G thin blind sterile single use cannula. The product is in ready-to-use form. The amount is specified on the product packaging. Products are sterilized by steam sterilization.

3.INTENDED PURPOSE

Corneal endothelial cells may be damaged, in case the space between the corneal endothelium and the anterior chamber is not maintained during surgery, causing corneal decompensation. Viscoelastic solution is injected to protect the corneal endothelium at the cavity in the anterior chamber. The main objective of viscoelastic is to create a protective barrier for the corneal endothelium or to ensure intraocular lenses to slide out of their cartridges during operation. Methylcellulose protects chamber depth, endothelial layer, other intraocular tissues during anterior segment surgery and this allows more efficient manipulation with significantly less trauma to the corneal endothelium and other ocular tissues.

3.1.Clinical Benefit

•Supports intraocular manipulations by creating volume and protecting the endothelium.

•Enables surgeons to easily perform surgical maneuvers during operations.

4.INDICATIONS

Methylcellulose solution is an ophthalmic surgical instrument used in surgical interventions regarding the anterior segment like cataract and intraocular lens implantation (IOL), routine cataract surgeries, intraocular lens implantation, corneal transplantation, and glaucoma surgeries.

4.1.Patient Information

The targeted patient group includes all age groups, all nationalities, and the whole population in need of ophthalmologic surgery without indicating any weight restrictions. The International Patient Identification Card is issued to the patient with instructions to keep the implant as a permanent record and to show it to any future ophthalmologist.

4.2.Patient Preparation

The product is used in eye surgery. The patient is prepared for the surgical procedure in line with the doctor's instructions. The product does not require any special preparation and training.

5.CONTRAINDICATIONS

The surgeon should use careful pre- and peri-operative evaluation and solid clinical judgment to decide the benefit/risk ratio before administering a methylcellulose injection to a patient with one or more of the following conditions. The evaluation of each individual case is at surgeon's discretion.

• In patients who are sensitive to any substance in the solution,

• In patients with very high intraocular pressure,

• Use for children, pregnant or lactating women is under the initiative of the ophthalmic surgeon.

6.SIDE EFFECTS

Temporary high intraocular pressure abnormalities may occur following the surgery. The patient should be warned that he/she should consult his/her doctor if intraocular pressure increases. Even though there are no cases reported up to the present, it may interact with concomitant medications and cause condensation. In case such a case is observed, physicians should consider this potential. Risk analysis conducted by the manufacturer's risk management team specified that the remaining risk is as low as reasonably possible.


7.USAGE

These are the instruments used during ophthalmic surgery. Should not be left in the eye after the operation is completed. The amount to be injected varies in each specific case and should be designated by the surgeon during the operation. Do not mix with any other substance before use.

•Take out the products from the box.

•Open the product packaging just before use.

•Open the sterile packaged products taken out of the box to a sterile area by following the sterile opening rules and the opening direction shown on the package. **(Picture A)**


Picture A: Symbol of Roll Pack Open Direction

Instructions to remove the syringe luer cap correctly:

1. Take the syringe from the box out of its package.
2. Hold the syringe luer as shown.
3. Carefully twist the luer cap.
4. Remove the luer cap.

Instructions to place the cannula correctly:

5. Take the cannula from the box out of its package.
6. Hold the syringe.
7. Place the cannula on the syringe.
8. Twist and secure.
9. After the air is purged out of the cannula tip, the solution is ready to be injected. Methylcellulose viscoelastic is completely cleaned/removed from the eye by irrigation and/or aspiration method at the end of the surgery.

Manufacturer: TEKNOMEK Medikal Malzemeleri Sanayi ve Ticaret Ltd. Şti. Karamehmet Mah. Avrupa Serbest Bölgesi 9. Sk. No:14 Ergene / Tekirdağ / Turkey Tel: +90 212 369 1263, Fax: +90 212 210 3432 Website: www.teknomek.net e-mail: info@teknomek.net Document No: IFU15M Rev.07/08.12.2023 (You can find the recent version and revision information on Teknomek's website.)

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Address: Kervangeçmez Cad. Hark Sok. Çağdaş Apt.

No:2/5 34381 Mecidiyeköy, Istanbul, Türkiye

Tel: +90 212 217 92 06

e-mail: info@atlantiktercume.com

14.SYMBOLS

The symbols given on the product box and in the instructions for use are listed below.

LOT	Manufacturer
25°C 15°C	Batch Code
30% 70%	Temperature Limit
	Humidity Limitation
	Do not re-use.
	Do not resterilize.
	Use-by Date
	Caution
	Keep away from sunlight
STERILE	Sterilized using steam. Single sterile barrier system.
	Keep dry.
	Consult instructions for use.
	Do not use if package is damaged and consult instructions for use.
MD	Medical device
UDI	Unique Device Identifier
	Does not contain natural rubber latex.

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1.PRODUKTBESCHREIBUNG

Produktname: Hydroxypropylmethylcellulose

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