445,5 mm

49,5 mm

200 mm

TEKNOMEK

TR

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OPHTHALMIC CANNULES AND NEEDLES INSTRUCTIONS FOR USE
Read the instructions for use cautiously!

1. PRODUCT INFORMATION
Product Name: Ophthalmic Cannulas and Needles
Product Model: Visit the catalog or website for the models
BUDI-DI: 868054899000N3
Shelf Life: Varies according to the cannula types. Stated on the product.
Class: IIa
Ophthalmic cannulas are sterile, single use, surgically invasive medical devices used mainly for irrigation, aspiration, and removal of fluids in eye operations. Ophthalmic needles are sharp-edged, sterile, single use, surgically invasive medical devices designed to be used for various purposes in eye surgeries with liquid substances.
The product comprises of stainless steel and hub. There are varieties containing polyimide and silicone tubing depending on the type.
Cannulas and needles are compatible for the syringe connection. Cannula gauge information varies depending on the hub color. **(Table 1: Gauge-Color Types)**

Table 1: Gauge-Color Types

COLOR	GAUGE	COLOR	GAUGE	COLOR	GAUGE
CREAM	19G	YELLOW	20G	GREEN	21G
BLUE	23G	PURPLE	24G	ORANGE	25G
BROWN	26G	GREY	27G	NATURAL	30G

2.PACKAGING and STERILIZATION
Products are sterilized by ethylene oxide sterilization. The product is offered sterile in blister or sterilization roll. Product information is stated on the product packaging.

3.INTENDED PURPOSE
Cannulas are used in anaesthesia, nucleus removal, irrigation and aspiration, capsule polishing, gaining access to the lens, cortex removal, hydro dissection, anterior chamber forming and drug administration in cataract, glaucoma, strabismus and refractive operations and vitreoretinal surgeries, aspiration and irrigation, injection, anterior chamber irrigating and cornea moistening.The type of cannula is chosen depending on the application purpose.

4.INDICATIONS
Anesthetize, creating an opening from the anterior capsule, administering viscoelastic and washing the cornea during the cataract surgery, separating the nucleus and cortex in cataract operation, irrigation and expansion of the lacrimal canal, forming the globe volume by infusion in vitreoretinal surgery, providing the outflow of in vitro fluid simultaneously while administering the decalin product, and operations of aspiration of substance from the subretinal areas are indicated.
4.1.Patient Information
The targeted patient group includes all age groups, all nationalities, and the whole population when needed for vitrectomy surgeries without indicating any weight restrictions.
4.2.Patient Preparation
The product is used in eye surgery. The patient is prepared for the surgical operation in line with the doctor's instructions. The product does not require any special preparation and training.

5.CONTRAINDICATIONS
Should not be used by patients known to be allergic to any of the product ingredients. The evaluation of each individual case is at surgeon's discretion.

6.SIDE EFFECTS
Risk analysis conducted by the manufacturer's risk management team specified that the remaining risk is as low as reasonably possible.

7.USAGE
•Take out the products from the box.
•Open the product packaging just before use.
•Open the sterile packaged products taken out of the box to a sterile area by following the sterile opening rules and the opening direction shown on the package.
(Picture A) (Picture B)

Picture A: Blister Packaging Opening Direction Symbol

Picture B: Roll Packaging Opening Direction Symbol

Picture C

•Install the product by transforming it into the system suitable for the operation to be performed. **(Picture C)**
•The product is ready to use

8.WARNINGS and PRECAUTIONS

It is for single use.
Do not re-sterilize. Reuse and/or resterilization may adversely affect device performance, which could seriously harm patient health and safety.
Do not use the product whose sterile package has been opened or damaged.
Do not use expired products.
The product is restricted for use by physicians specialized in surgery trained in ophthalmic surgery.
It must be considered that the tip of the product may be sharp, care should be taken while using.
Stop the application if any side effects occur. In case of any adverse reaction, notify the manufacturer.
Conduct the application according to the instructions written in the instructions for use. Otherwise, the manufacturer is not responsible.
Use for children, pregnant or lactating women is under the initiative of the ophthalmic surgeon.
Pay attention to the features of use as per the type of product you buy.

9.WASTE DISPOSAL
Used products must be disposed of properly to avoid possible risks of infection. Disposal of the products after use should be performed by authorized persons in line with hospital procedures. Manufacturer is not responsible for waste disposal procedures.

10.STORAGE AND TRANSFER CONDITIONS
Store at temperatures between 10°C (50°F or 283.15°K) and 35°C (95°F or 308.15°K). Store in a dry place. Protect against light and impact. If it is wanted to be transferred, it can be transported in the boxes belonging to the product.

11.WARRANTY AND LIMITATION OF LIABILITY
Teknomek Medikal is not responsible for any injuries or damages that may arise from non-compliance with the recommendations and suggestions in the instructions for use and that are not suitable for the patient in the specified criteria. The risk of contamination which may arise from re-sterilizing the product by ignoring its single use quality is not the responsibility of the manufacturer. Teknomek Medikal is not responsible for the injuries or damages that the patient may suffer if it is not used by the doctor in line with the instructions for use.

12.SAFETY and CLINICAL PERFORMANCE
The product does not contain tissue, cell or derivatives of human or animal. There are no accessories and software for the device. There are no medicinal substances or biological materials included in the device as an integral part of the device. No materials that contain or consist of endocrine disruptors as an integral part of the device or that may cause sensitization or an allergic reaction to the patient or user are included in the device. The products contain trace amounts of Nickel, which is one of the CMR substances. As a result of the risk management report conducted by the manufacturer's risk management team, the toxicological risks that the product may transfer to the patient in clinical use have been evaluated specific to contact time and the tissues where the contact is made, and it shows that the product is not cytotoxic, ocular irritant, skin sensitive and acute systemic toxic. In line with this result, the product is biocompatible. Please contact the manufacturer for detailed safety and performance information.

13.RETURN AND EXCHANGE POLICY
Please contact the manufacturer or your local distributor for the return or exchange of the product.

14.SYMBOLS
The symbols given on the product box and in the instructions for use are listed below.

	Manufacturer
LOT	Batch Code
-35°C 10°C	Temperature Limit
70% 20%	Humidity Limitation
	Do not re-use.
	Do not resterilize.
	Use-by Date
	Caution
	Keep away from sunlight
STERILE EO	Sterilized using ethylene oxide. Single sterile barrier system.
	Keep dry.
	Consult instructions for use.
	Do not use if package is damaged and consult instructions for use.
MD	Medical device
UDI	Unique Device Identifier
	Does not contain natural rubber latex.

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