49,5 mm

594 mm

200 mm

TEKNOMEK

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SILICONE OIL INSTRUCTIONS FOR USE

Read the instructions for use cautiously!

EN

1.PRODUCT INFORMATION

Product Name: Silicone Oil
Product Model: See the catalog or web-site for the models
BUDI-DI : 868054899021NB
Size: 10 ml (in glass vial or glass syringe)
Shelf Life: Shelf life is 2 (two) years. Stated on the product
Class: IIb

Technical	
Boiling Point:	>205 C
Melting Point:	<-40 C
Specific Weight:	0.971
Smell:	Odorless
Water Solubility	Insoluble (<100ppb)


1.1 Chemical Properties

Silicone Oil	1000 cSt
Density	0.961 g/cm³- 0.970 g /cm³
Viscosity	1000 cSt
Endotoxicity	<0.2EU/ml
Retractive Index	1.403-1.404 at 25 ˚C
Spectral Transmittance	Given graphically between 300 nm and 1100 nm.

Silicone Oil	1300 cSt
Density	0.961 g/cm³- 0.970 g /cm³
Viscosity	1300 cSt
Endotoxicity	<0.2EU/ml
Retractive Index	1.403-1.404 at 25 ˚C
Spectral Transmittance	Presented graphically between 300 nm and 1100 nm.

Silicone Oil	5000 cSt
Density	0.965 g/cm³- 0.973 g /cm³
Viscosity	5000 cSt
Endotoxicity	<0.2EU/ml
Retractive Index	1.403-1.404 at 25 ˚C
Spectral Transmittance	Presented graphically between 300 nm and 1100 nm.

Spectral Transmittance Curve (1000 cSt)

Spectral Transmittance Curve (1300 cSt)

Spectral Transmittance Curve (5000 cSt)

Silicone oil is a transparent, odorless, colorless, sterile, single use, surgically invasive, ultrapure medical liquid. It is a purified, sterile vitreoretinal, fractionated apyrogenic polydimethylsiloxane silicone oil for surgical purposes. This product does not contain phthalates and latex. Silicone oil is a very pure substance. Silicone oil is used as vitreous tamponade in vitrectomy surgery. For an easier application, silicone oil can be supplied with syringe, cannula, injection system.

2.PACKAGING and STERILIZATION

Products are sterilized by steam sterilization. The products are offered in glass vial sealed with a polymer stopper and metal cap or in a syringe as a 10 ml. Depending on the disposable product code, the silicone injection cannula is supplied with syringe and small glass vial form. The syringe form is ready to use. The amount and viscosity are stated on the product package.

3.INTENDED USE

Silicone oil is used as a substitute for the vitreous body in vitrectomy surgery. The product is preferred for a long-term ocular endotamponade in selected areas of complicated or uncomplicated recurrent retinal detachments.

3.1.Clinical Benefit

Helps in the surgical treatment of retinal detachment as a result of giant tear retinopathy, proliferative vitreoretinopathy, fractional diabetic retinopathy and traumatic retinopathy and other viral infections affecting the retina

4.INDICATIONS

It is a substance used as retinal tamponade (up to 29 days) postop with all RD (giant tear retinopathy, proliferative vitreoretinopathy, fractional diabetic retinopathy and traumatic retinopathy).

4.1.Patient Information

The targeted patient group includes all age groups, all nationalities, and the whole population acquiring endotamponade following the vitrectomy surgery without specifying any weight restrictions. The international patient identification card is given to the patient with instructions to keep the implant as a permanent record and to show it to any ophthalmologist in future.

4.2.Patient Preparation

The product is used in eye surgery. The patient is prepared for the surgical procedure in line with the doctor's instructions. The product does not require any special preparation and training.

5.CONTRAINDICATIONS

Surgeon should use careful pre- and peri-operative evaluation and solid clinical judgment to decide the benefit/risk ratio before injecting silicone oil to a patient with one or more of the following conditions. The evaluation of each individual case is at surgeon's discretion.

- In patients known to be allergic to any of the product ingredients,
- In patients with known silicone allergy,
- In pseudophakic patients with a silicone intraocular lens (IOL).

6.SIDE EFFECTS

There are also risks in surgeries performed using silicone oil for the superior segment as with all surgical interventions. Usage below the required amount poses the risk of ineffective and low tamponade. Redundant use may cause problems such as passing to anterior chamber and corneal problems, cataracts, and secondary glaucoma.

Ocular hypertonia, ocular hypotonia, choroidal detachment, corneal epithelial edema, cystoid macular edema, infection, endophthalmitis, glaucoma, keratopathy, transfer of silicone oil to the anterior chamber, anti silicone antibodies and opacification have been reported among the silicone-related side effects. Other side effects include optic nerve atrophy, rubeosis iridis, transient intraocular pressure increase, macular pucker, fibrosis and rupture, which are observed less. Passing of silicone oil into the anterior chamber may result in corneal changes including opacity, descemet's folds, or corneal decompensation. If iridectomy is not performed at 6 o'clock position in aphakic eyes, pupillary block and angle closure may occur after silicone oil infusion. Risk analysis conducted by the manufacturer risk management team specified that the remaining risk is as low as reasonably possible.

7.USAGE

The amount to be injected will vary in each specific case and must be determined by the surgeon during the operation. It must be checked whether there is a problem in the correct papillary circulation and feeding of the optic nerve during the infusion because it is impossible to determine the correct intraocular pressure in practice. In case suspected for such a situation, the infusion must be stopped immediately. Do not mix with any other substance before use.

- Take out the products from the box.
- Open the product packaging just before use.
- Open the sterile packaged products taken out of the box to a sterile area by following the sterile opening rules and the opening direction shown on the package. **(Picture A)**
- Be attention not to form air bubbles while drawing into the syringe from the glass vial.

Picture A: Roll Pack Open Direction

7.1. Correct Instructions for Use for Glass Vial:

1. Take the glass vial from the box out of its package.
2. Remove the cap of the glass vial.
3. Take the cannula from the box out of its package.
4. Hold the syringe.
5. Place the cannula on the syringe (which is suitable for delivery)
6. Twist and secure.
7. Perforate the polymer plug with the syringe. Pull the plunger rod to draw the product into the syringe.
8. Transfer the product from the vial to the syringe

1. Take the syringe out of its package with its plunger rod.

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